

# Invacare® IVC™ 9000 XDT Wheelchair

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### SEAT-TO-FLOOR HEIGHT

These selections are based on 24" wheels and 8" casters.

### TRANSACTION TYPE:

Order  Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

## FRAME AND SEAT SELECTIONS

### FRAME TYPE

See seat-to-floor height selections.

9XDT 9000 XDT Heavy-Duty (K0007) ..... \$1,544.38 K0007

### FRAME WIDTH AND DEPTH

- WD66 16" x 16" ..... No Charge
- WD67 16" x 17" ..... No Charge A9900
- WD68 16" x 18" ..... No Charge A9900
- WD60 16" x 20" ..... \$126.38 E2203
- WD76 17" x 16" ..... No Charge A9900
- WD77 17" x 17" ..... No Charge A9900
- WD78 17" x 18" ..... No Charge A9900
- WD70 17" x 20" ..... \$126.38 E2203
- WD86 18" x 16" ..... No Charge
- WD87 18" x 17" ..... No Charge A9900
- WD88 18" x 18" ..... No Charge A9900
- WD89 18" x 19" ..... \$126.38
- WD80 18" x 20" ..... \$126.38 E2203
- WD96 19" x 16" ..... No Charge A9900
- WD97 19" x 17" ..... No Charge A9900
- WD98 19" x 18" ..... No Charge A9900
- WD99 19" x 19" ..... \$126.38
- WD90 19" x 20" ..... \$126.38 E2203
- WD06 20" x 16" ..... No Charge E2201
- WD07 20" x 17" ..... No Charge E2201
- WD08 20" x 18" ..... No Charge E2201
- WD09 20" x 19" ..... \$126.38 E2201
- WD202 20" x 20" ..... \$126.38 E2201 E2203
- WD26 22" x 16" ..... No Charge E2201
- WD27 22" x 17" ..... No Charge E2201
- WD28 22" x 18" ..... No Charge E2201
- WD29 22" x 19" ..... \$126.38 E2201
- WD222 22" x 20" ..... \$126.38 E2201 E2203

### SEAT-TO-FLOOR HEIGHT

These selections are based on 24" wheels and 8" casters.

- A Adult (19 1/2" STF) ..... Standard
- H Hemi (17 1/2" STF) ..... No Charge
- T Tall (21 1/2" STF) ..... No Charge K0056

## ARM AND BACK OPTIONS

### ARM TYPE

These selections come standard with removable flipback arm styles.

- 18 Full Length Fixed Ht Conv (10") ..... No Charge A9900
- 19 Full Length Adj Ht Conv (10"-14") ..... \$134.95 E0973
- 28 Desk Length Fixed Ht Conv (10") ..... No Charge A9900
- 29 Desk Length Adj Ht Conv (10"-14") ..... \$134.95 E0973
- 19JY Full Length Adj Ht Conv (8"-10") ..... \$139.23 E0973
- 29JY Desk Length Adj Ht Conv (8"-10") ..... \$139.23 E0973
- 78 Full Length Fixed Ht SS (10") ..... No Charge
- 79 Full Length Adj Ht SS (10"-14") ..... \$134.95 E0973
- 88 Desk Length Fixed Ht SS (10") ..... No Charge
- 89 Desk Length Adj Ht SS (10"-14") ..... \$134.95 E0973

### BACK TYPE

- ADJ Adjustable Back with 10° Bend 15"-19" ..... No Charge
- S1418 Straight Back with Push Handles (14"-18")... No Charge K0108

## COLOR OPTIONS

### FRAME FINISH

- 17P Cobalt Blue ..... No Charge
- 61P Electric Red ..... No Charge
- 24P Wet Black ..... No Charge
- 34P Satin (flat) Black ..... Standard
- 71P Silver Vein ..... No Charge

### UPHOLSTERY COLORS

- U240 Black Nylon Upholstery ..... Standard
- U245 Black Sail Cloth Adj Tension Seat ..... No Charge  
**U245 sail cloth adjustable tension seat upholstery available with U240 black nylon back upholstery only.**
- U67 Black Vinyl ..... No Charge  
**Vinyl Upholstery must choose WRAP or SOUPH.**

### NAUGAHYDE

- WRAP Wrap Around Vinyl Upholstery ..... \$82.47 A9900
- SOUPH Screwed On Back Upholstery ..... \$95.32 A9900

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

# Invacare® IVC™ 9000 XDT Wheelchair

## CASTERS

### CASTER OPTIONS

- 1228 8" x 1.75" Semi Pneu Casters Light Grey ..... \$59.98 A9900
- AT903 8" x 1.25" Comp Urethane ..... Standard
- 1299 8" x 1" Solid Rubber Casters ..... No Charge

## WHEELS

### WHEEL OPTIONS

- 5244BK 24" HD Comp Pneu W/FF Inserts ..... \$319.16 E2213
- U2244BK 24" Comp Urethane Treaded Tires ..... \$128.52
- U2222C 24" Comp Urethane Corded Tires ..... Standard

### WHEEL OPTIONS (ONE ARM DRIVE)

These selections are not available with amputee attachment.

- 9456L One Arm Drive with U2222C LH ..... \$649.03 E0958
- 9456R One Arm Drive with U2222C RH ..... \$649.03 E0958

## HANDRIM OPTIONS

- AHDRM Aluminum Hndrms ..... No Charge  
**Not available with standard U2222C rear wheels.**
- 1512 Plastic Coated Aluminum Hndrms ..... \$91.04 A9900  
**Not available with 5244BK or U2244BK wheel options.**
- 1224-8 8 Alum Oblique Proj Hndrms, Plst Coated ..... \$269.89 A9900
- 1225-8 8 Alum Vert Proj Hndrms, Plst Coated ..... \$269.89 A9900

## WHEEL LOCKS

- 9690B Pull-To-Lock Pivot Link Whl Locks (Std Hemi) .. No Charge
- 9691B Push-To-Lock Pivot Link Whl Locks (Std Adult) . No Charge

## FRONT RIGGING OPTIONS

### FRONT RIGGING

- T93HA Hemi S/A Frst Alum Ftplt with Heel Loops ..... \$99.13 E0951
- T93HC Hemi S/A Frst Comp Ftplt with Heel Loops ..... \$99.13 E0951
- T94HA Hemi Elev Legrests Aluminum Ftplt ..... \$191.86 K0195 E0990
- T94HC Hemi S/A Elev Legrests Comp Ftplt ..... \$191.86 K0195 E0990
- AHL4A Hemi Smartleg® Art Legrests 16" and Above ... \$341.09 K0053
- AHL4C Hemi Smartleg® Art Legrests 16" and Above ... \$319.77 K0053
- T93HAB Hemi Bolt-Thru Footrests, Aluminum Footplates ..... \$95.93 A9900
- T94HAB Hemi Bolt-Thru Legrests, Aluminum Footplates ..... \$186.53 A9900

## FOOTPLATE OPTIONS

- 1350 Extra Large Aluminum Footplates ..... \$84.21 K0041
- AT5543 Adjustable Angle Flip-Up Footplates ..... \$71.42 K0040

## FOOTREST/LEGREST ACCESSORIES

- 1820 Non-Removable Device ..... \$12.79 A9900
- 1337 Calf Strap ..... \$38.37 K0038
- 1346 Heel Loops ..... \$41.57 E0951  
**Standard on Swingaway footrests.**

## POSITIONING STRAP OPTIONS

### SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap ..... \$53.55 A9900 E0978
- 1321BK Hook and Loop Seat Positioning strap ..... \$42.84 E0978

## ACCESSORIES

### BACK ACCESSORIES

- L8601 Omit Back Upholstery ..... Less \$0.00
- 1160 Carrying Pocket ..... \$57.83 A9270
- STC Stencil Back Upholstery ..... \$31.06 A9270  
**Only available with vinyl upholstery. 3/4" block letters.**
- RB Dual Back Liner ..... \$49.00 A9270

### SEAT ACCESSORIES

- L8651 Omit Seat Upholstery ..... Less \$0.00
- 9120 Solid Seat Insert ..... \$104.96 A9900
- RS Dual Seat Liner ..... \$52.48 A9270

### ARM ACCESSORIES

- NRAS Non-Removable Arm with Screws ..... \$10.71 A9900
- 1300 Utility Tray ..... \$168.15 E0950  
**Not available with 24" rear wheel set in Hemi position.**

### MISCELLANEOUS ACCESSORIES

- 9357L Amputee Attachment ..... \$100.78 E0959  
**Not available with one arm drive.**
- 1305L Crutch and Cane Carrier ..... \$46.05  
**Not available with one arm drive.**
- 2226 Front Anti-Tipper ..... \$88.89  
**Only available with H93M front riggings.**

### STOCK OPTIONS

- 1100 Tool Kit ..... No Charge
- 1497 Oxygen Holder Package Assembly ..... \$161.72 E2208  
**Only available on Adult or Hemi seat-to-floor.**
- 1496 IV/Oxygen Holder Package Assembly ..... \$169.22 K0108  
**Only available on Adult or Hemi seat-to-floor.**
- 1360A Adjustable Rear Anti-Tipper ..... \$87.82 E0971
- A1576 Wheel Lock Extensions ..... \$41.77 E0961  
**Elastic cord wheel lock extension. Seat-to-floors come standard with non-cord wheel lock extension.**


All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

**Invacare® IVC™ 9000 XDT Wheelchair**

<b>Seat-to-Floor 8" Casters</b>			
<b>Seat-to-Floor Height</b>	<b>Front Caster Mounting Position</b>	<b>Rear Wheel Size</b>	<b>Rear Wheel Mounting Position</b>
17 1/2"	Top	24"	Top
19 1/2"	Top w/2-inch extension	24"	Middle
19 1/2"	Bottom	24"	Middle
21 1/2"	Bottom w/2-inch extension	24"	Bottom



**Weight Limitation 350 lbs.**  
 Invacare recommends all wheelchairs to be fitted to the consumer under the guidance of a healthcare professional knowledgeable of the particular limitations of the product.

Comments

**All prices subject to change without notice.**  
**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*

**Invacare Corporation**

[www.invacare.com](http://www.invacare.com)

**USA**

One Invacare Way  
 Elyria, OH USA  
 44035-4190  
 800-333-6900

©2020 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.



Revised January 1, 2020