


Frame Assembly, Front and Rear



Frame Assembly, Front and Rear

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1099721-NLA	NLA - Strut, Housing and Hardware Includes: 4 ea. Strut Housing, 4 ea. O-Ring, Rubber, 4 ea. Snap Ring	4	4
2		1099717-NLA	NLA - Axle Assembly (Rear) Includes: 1 ea. Axle, Rear 1 ea. Bolt, Square Head (M12 x 1.25 x 68) 1 ea. Spacer (M12 x 20 x 4.5) 1 ea. Locknut (M12 x 1.25 x 19) 2 ea. Shroud Cover, Rear Axle 2 ea. Screw, Phillips Pan Head (M4 x 12)	1	1
3		1099724-NLA	NLA - Hardware, Connecting Pin Includes: 4 ea. Connecting Pin, Frame 4 ea. Screw, Phillips Self Tap (M2 x 10)	4	4
4		1099722-NLA	NLA - Hardware, Adjustable Stop Rear Frame Includes: 1 ea. Washer, Spring (M8 x 13) 1 ea. Nut (M8 x 1.25 x 13) 1 ea. Screw, Set (M8 x 1.25 x 35) 1 ea. Fastener, Hook Adhesive Back	1	1
5		1099716-NLA	NLA - Axle Assembly (Front) Includes: 1 ea. Axle, Front 1 ea. Bolt, Square Head (M12 x 1.25 x 68) 1 ea. Spacer (12 x 20 x 4.5) 1 ea. Locknut (M12 x 1.25 x 19) 1 ea. Shroud Cover, Front Axle - Right 1 ea. Shroud Cover, Front Axle - Left 2 ea. Screw, Phillips Pan Head (M4 x 15)	1	1