





# 16" Seat and Arms



## 16" Seat and Arms

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1117008-NLA	NLA - Seat Assembly, Grey (16" W)	1	1
2	C,D	1116977-NLA	NLA - Kit, Arm Assembly with Seat Pivot (16" Wide Seat)	1	1
3			NLA - Pivot, Seat	1	1
4			Washer, Steel (5/16 x 1-1/4 x 1/16")	1	4
5			NLA - Screw, Hex Head w/Patch (1/4-20 x 1-1/8")	1	4
6			NLA - Bushing (1/4 x 3/4 x 1/2")	1	6
7			NLA - Arm Assembly, Flip-Up, Black - Right	1	1
8			Knob, Adjustment, Black (5/16-18 x 3/4")	1	3
9			NLA - Arm Assembly, Flip-Up, Black - Left	1	1
10			NLA - Bracket, Arm Support	1	1
11			Washer (1/4 x 1/2 x 1/16")	1	6
12			NLA - Screw, Hex Head w/Patch (1/4-20 x 1-3/8")	1	6
	A,B	1070553-NLA	NLA - Screw, Phillips Pan Head (#10-32 x 5/8")	1	4
13	E	1114228-NLA	NLA - Arm Assembly, Flip-Up, Black - Left	1	1
13		1137221-NLA	NLA - Arm Assembly, Flip-Up, Black - Left	1	1
14		1116976-NLA	NLA - Arm Pad Assembly, Desk Length with Hardware - Left	1	1
15		1114910	Knob, Adjustment, Black (5/16-18 x 3/4")	1	1
16	E	1116906-NLA	NLA - Tube, Lower Arm - Left	1	1
16		1113647	Tube, Arm, Lower Adjustable Width, Black - Left	1	1
17	E	1114227-NLA	NLA - Arm Assembly, Flip-Up, Black - Right	1	1
17		1137220-NLA	NLA - Arm Assembly, Flip-Up, Black - Right	1	1
18		1116975-NLA	NLA - Arm Pad Assembly, Desk Length with Hardware - Right	1	1
19	E	1116905-NLA	NLA - Tube, Lower Arm - Right	1	1
19		1113646	Tube, Arm, Lower Adjustable Width, Black - Right	1	1

NOTE: A - Not Shown.  
B - Screws to attach Arm Pad.  
C - Items 7 & 9, 13 & 17 Arm Assembly Includes: Items 14 & 18 (Arm Pads) and item 15 (Adjustment Knobs).  
D - Includes: Items 2-12.  
E - Due to style change, Arm Rests need to be replaced in pairs. Order 1 ea. of 1137221 Left and 1137220 Right