


Rear Wheels - Mountain Bike Wheel - 24" Wheels

NOTE: 24" Mountain Bike Wheel Assembly uses a 22" Handrim.



Rear Wheels - Mountain Bike Wheel - 24" Wheels

NOTE: 24" Mountain Bike Wheel Assembly uses a 22" Handrim.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1562665	Rear Wheel (24"), Mountain Bike, Aluminum Handrim, Black, MTB Tire	1	2
2		1529327	Rear Wheel (24"), Mountain Bike, No Handrim, No Tire	1	2
3		1571073	Spoke Reflector (Pair)	2	1
4		1579085	Kit, Spoke, Mountain Bike Wheel (10 Pcs.)	10	1