



Headrest Extension - Non-Recliner



16" Wide Non-Reclining Chairs



14" Wide Non-Recliner Chairs

Headrest Extension - Non-Recliner

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	B,C	88801066X000-NLA	NLA - Upholstery, Headrest	1	1
1	B,D	8881054615-NLA	NLA - Upholstery, Headrest	1	1
2		1027095	Plug Button, Black (7/8")	1	2
3		1076283	Screw, Phillips Oval Head w/ Washer (#10-16 x 7/8")	1	6
3		1041541	Package, Screw Phillips Csunk Oval Head w/ Washer (#10-16 x 7/8") (12 ea.)	12	1
4	A,C	1011725-NLA	NLA - Headrest Extension Tube Assembly	1	2
4	D	01058X000-NLA	NLA - Tube, Headrest Extension	1	2
4	A,C	1054610-NLA	NLA - Headrest Extension Tube w/o Snap Fasteners Assembly	1	2
4	D	1054604-NLA	NLA - Tube, Headrest Extension	1	2
5	C	48811X004-NLA	NLA - Insert, Tube	1	2
5	A,D	48811X003	Tube, Insert	1	2
6		47447X009-NLA	NLA - Screw, Phillips (#8 x 3/8")	1	6
7		48127X002-NLA	NLA - Screw, Phillips Pan Head (#10-32 x 5/8")	1	1
8		04546X005-NLA	NLA - Bracket, Chrome	1	2
9	B,C	8881054614-NLA	NLA - Headrest Extension Assembly	1	1
9	B,D	1054613	Headrest Extension Assembly	1	1

NOTE: A - For 14" wide chairs built after 6-14-94, the insert tube is now part of the extension tube (item 4) and cannot be purchased separately.
B - Specify upholstery color, see page titled "Packaging and Upholstery Colors"
C - 14" Wide Chairs
D - 16" Wide Chairs