


Seat Pan & Extensions

Basic Recline 4XP

NOTE: To achieve an 18" seat depth, order a 16" deep seat pan and 2" extension kit.

NOTE: To achieve an 20" seat depth, order a 16" deep seat pan and 4" extension kit.



Item #1 Seat Pan

Depth	16" Wide	18" Wide
16" Deep	0221419	0221519

Seat Pan & Extensions

Basic Recline 4XP

NOTE: To achieve an 18" seat depth, order a 16" deep seat pan and 2" extension kit.

NOTE: To achieve an 20" seat depth, order a 16" deep seat pan and 4" extension kit.

Item #2 Seat Pan Extensions

Depth	16" Wide	18" Wide
2" Deep	0123019-NLA	0123219-NLA
4" Deep	0221119-NLA	0221219-NLA

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Seat Pan	1	1
2		See Chart	#2 Seat Pan Extensions	1	1
3		1072695-NLA	NLA - Screw, Phillips Pan Head Tap (#6 x 1/2")	1	4