

# Rollator (Four Wheeled) - Model 65650 and 65650R



## Rollator (Four Wheeled) - Model 65650 and 65650R

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1142628	Grip, Hand (Pair)	1	1
2	B,F,I	1142629	Brake Handle w/ Hardware Assembly	1	2
3	C,G,I	1142630	Brake Assembly	1	2
4		1142632	Back Support Assembly, Curved	1	1
5		1142633	Seat w/ Hardware	1	1
6		1142638	Basket	1	1
7		1142639	Connection Linkage Assembly	1	2
8		1142637	Package, Strap	10	1
9	D	1142636	Kit, Wheel w/ Hardware - Rear	1	2
10	H	1142634	Fork Assembly, Front	1	2
11		1142635	Kit, Wheel w/ Hardware - Front	1	2
12	E	1144495-NLA	NLA - Kit, Adjustable Knob (Pair)	2	1
12		1169007	Kit, Adjustable Knob (Pair)	2	1

NOTE: A - Includes (2) Grips.  
 B - Includes Handle with attaching Hardware.  
 C - Includes Handle with Cable, Fork, Brake Pad and Hardware.  
 D - Includes Rear Wheel with Hardware.  
 E - Part is no longer available. See #1169007 for replacement.  
 F - Part #1142629 is used for Model 65650 and 65650R manufactured before 11/01/07 only.  
 G - Must purchase as a pair for all Rollators manufactured before 4/16/2008.  
 H - Replacement wheel has different amount of Spokes that does not impact performance. It is recommended to replace Fork Assembly in pairs to maintain appearance.  
 I - New Handgrips must be ordered as existing Handgrips must be cut off to install new Brake Assembly.